



# Contec Sterile Alcohol Low Endotoxin

**70% IPA or denatured ethanol blended with WFI or purified water**

Contec alcohols are a blend of 70% v/v Isopropanol or denatured ethanol with 30% water for injection or purified water.

The alcohol blend is 0.2 micron filtered, filled and bagged in a Grade C (ISO Class 7) cleanroom. This clean manufacture coupled with water for injection or purified water means the alcohol blend is guaranteed to have an endotoxin level of less than 0.25EU / ml.

Contec's ethanol is only denatured with IPA creating a very low residue product. Contec sterile alcohol is provided sterile by gamma irradiated using a validated process at no less than 25 kGy.

Supplied as 0.5L or 1L trigger sprays fitted with a protected system, which ensures sterility throughout use, or 5L capped container for larger areas.



FEATURES	BENEFITS
Guaranteed endotoxin levels of <0.25EU/ml	• Suitable for use in product contact areas
GMP manufactured under Grade A laminar flow in a Grade C cleanroom	• Ensures the alcohol, container and packaging are free from contamination and particulates
No denaturants other than IPA	• Completely residue free
Trigger sprayer and "bag in bottle" protected system	• Prevents bottle contents from becoming contaminated during use • Bottles can be completely emptied eliminating wastage so no need to discard unused product during shelf life
Trigger sprayer can be set to jet or spray	• Large droplet size reduces the risk of inhalation and provides good surface coverage
Double/triple bagged packed in linear tear packaging	• Each bag is easy to open even when wearing gloves • Facilitates transfer disinfection into high grade areas
Sterile via gamma irradiation	• Suitable for use in Grade A and B cleanrooms

For more information or to request a sample, e-mail [infoeu@contecinc.com](mailto:infoeu@contecinc.com) or phone + 33 (0) 2 97 43 76 98.

Copyright© 2019 Contec, Inc. All rights reserved. INT124 061119

Contec  
Vannes Cedex  
France  
Tel: +33(0) 2 97 43 76 98

Contec Inc  
Spartanburg, SC  
United States  
Tel: +1 (864) 503-8333

Contec Cleanroom Technology (Suzhou) Co. Ltd  
Suzhou  
China  
Tel: +86-512-6274 4050



## Product Specification

### Trigger Spray

70% IPA or denatured ethanol with 30% WFI  
 Ready to use with a 3 year shelf life  
 0.2 micron filtered  
 Guaranteed endotoxin level of < 0.25 EU/ml  
 Tripped bagged in polyethylene linear tear bags  
 Tamper evident triggers  
 Sterile - irradiated at no less than 25 kGy

### Capped

70% IPA or denatured ethanol with 30% purified water  
 Ready to use with a 3 year shelf life  
 0.2 micron filtered  
 Guaranteed endotoxin level of < 0.25 EU/ml  
 Double bagged in polyethylene linear tear bags  
 Tamper evident caps  
 Sterile - irradiated at no less than 25 kGy

Thermoplastic polymers used which show good material compatibility with hydrogen peroxide and HPV

Certificates of analysis, sterility, irradiation and quantitative endotoxin levels, less than 0.25 EU/ml provided with every batch.

Full supporting documentation including efficacy information (EN1276 and EN1650) is available on request

## Product Application

Ready to use so can be applied directly to the surface ensuring even coverage. Physically remove contaminants from the surface using an appropriate cleanroom wipe or mop and a recommended wiping technique for optimum contamination control.

Part No.	Description	Size	Packaging
SBT0570IW	Contec Sterile IPA 70% IPA with water for injection, trigger spray	0.5L	8/case
SBT170IW	Contec Sterile IPA 70% IPA with water for injection, trigger spray	1L	6/case
SBC570I	Contec Sterile IPA 70% IPA with purified water, capped	5L	2/case
SBT0570DEW	Contec Sterile Denatured Ethanol 70% DE with water for injection, trigger spray	0.5L	8/case
SBT170DEW	Contec Sterile Denatured Ethanol 70% DE with water for injection, trigger spray	1L	6/case
SBC570DE	Contec Sterile Denatured Ethanol 70% DE with purified water, capped	5L	2/case

Contec offers a range of innovative wipes and mopping systems designed for efficient and effective cleaning in controlled environments. See [www.contecinc.com](http://www.contecinc.com) or contact your local sales representative for advice and information about the most suitable products and cleaning techniques for your facility.

*Use biocides safely. Always read the label and product information before use.*